

Product datasheet for **GA116460**

Human **FAM123C (AMER3)** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	AMER3 CRISPRa kit - CRISPR gene activation of human APC membrane recruitment protein 3
Symbol:	AMER3
UniProt ID:	Q8N944
Synonyms:	FAM123C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)