

Product datasheet for **GA116463**

Human DIPK2A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DIPK2A CRISPRa kit - CRISPR gene activation of human divergent protein kinase domain 2A
Symbol:	DIPK2A
Components:	GA116463G1, DIPK2A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116463G2, DIPK2A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116463G3, DIPK2A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8NDZ4
Synonyms:	DIA1; GoPro49; HASF
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)