

Product datasheet for **GA116465**

Human KIAA2018 (USF3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	USF3 CRISPRa kit - CRISPR gene activation of human upstream transcription factor family member 3
Symbol:	USF3
Components:	GA116465G1, KIAA2018 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116465G2, KIAA2018 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116465G3, KIAA2018 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q68DE3</u>
Synonyms:	KIAA2018
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

