

Product datasheet for **GA116514**

Human Connexin 40.1 (GJD4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	GJD4 CRISPRa kit - CRISPR gene activation of human gap junction protein delta 4
Symbol:	GJD4
Components:	GA116514G1, Connexin 40.1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116514G2, Connexin 40.1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116514G3, Connexin 40.1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q96KN9</u>
Synonyms:	CX40.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)