

Product datasheet for **GA116551**

Human MS4A15 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MS4A15 CRISPRa kit - CRISPR gene activation of human membrane spanning 4-domains A15
Symbol:	MS4A15
Components:	GA116551G1, MS4A15 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116551G2, MS4A15 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116551G3, MS4A15 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8N5U1
Synonyms:	FLJ34527; MGC35295
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)