

Product datasheet for **GA116566**

Human DOK6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DOK6 CRISPRa kit - CRISPR gene activation of human docking protein 6
Symbol:	DOK6
Components:	GA116566G1, DOK6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116566G2, DOK6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116566G3, DOK6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PKX4</u>
Synonyms:	DOK5L; HsT3226
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)