

Product datasheet for **GA116569**

Human HEPACAM activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HEPACAM CRISPRa kit - CRISPR gene activation of human hepatic and glial cell adhesion molecule
Symbol:	HEPACAM
UniProt ID:	<u>Q14CZ8</u>
Synonyms:	GlialCAM; MLC2A; MLC2B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)