

Product datasheet for **GA116571**

Human TIGD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TIGD3 CRISPRa kit - CRISPR gene activation of human tigger transposable element derived 3
Symbol:	TIGD3
UniProt ID:	<u>Q6B0B8</u>
Synonyms:	tigger transposable element derived 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)