

Product datasheet for **GA116577**

Human ZNF438 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF438 CRISPRa kit - CRISPR gene activation of human zinc finger protein 438
Symbol:	ZNF438
UniProt ID:	Q7Z4V0
Synonyms:	bA330O11.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)