

Product datasheet for **GA116585**

Human REEP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	REEP3 CRISPRa kit - CRISPR gene activation of human receptor accessory protein 3
Symbol:	REEP3
Components:	GA116585G1, REEP3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116585G2, REEP3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116585G3, REEP3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NUK4</u>
Synonyms:	C10orf74; Yip2b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)