


Product datasheet for **GA116593**

Human HNRNPUL2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HNRNPUL2 CRISPRa kit - CRISPR gene activation of human heterogeneous nuclear ribonucleoprotein U like 2
Symbol:	HNRNPUL2
Components:	GA116593G1, HNRNPUL2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116593G2, HNRNPUL2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116593G3, HNRNPUL2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q1KMD3</u>
Synonyms:	HNRNPUL2; SAF-A2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)