

Product datasheet for **GA116612**

Human GPR115 (ADGRF4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ADGRF4 CRISPRa kit - CRISPR gene activation of human adhesion G protein-coupled receptor F4
Symbol:	ADGRF4
UniProt ID:	Q8IZF3
Synonyms:	GPR115; PGR18
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)