

Product datasheet for **GA116614**

Human TDRD6 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TDRD6 CRISPRa kit - CRISPR gene activation of human tudor domain containing 6
Symbol:	TDRD6
Components:	GA116614G1, TDRD6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116614G2, TDRD6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116614G3, TDRD6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O60522</u>
Synonyms:	bA446F17.4; CT41.2; NY-CO-45; SPATA36; TDR2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)