

Product datasheet for GA116617

Human C6ORF154 (LRRC73) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC73 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 73
Symbol:	LRRC73
UniProt ID:	Q5JTD7
Synonyms:	C6orf154
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)