


## Product datasheet for **GA116645**

### Human PPP1R35 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PPP1R35 CRISPRa kit - CRISPR gene activation of human protein phosphatase 1 regulatory subunit 35
Symbol:	PPP1R35
Components:	GA116645G1, PPP1R35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116645G2, PPP1R35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116645G3, PPP1R35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TAP8</u>
Synonyms:	C7orf47
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)