

## Product datasheet for **GA116655**

### Human TMED4 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TMED4 CRISPRa kit - CRISPR gene activation of human transmembrane p24 trafficking protein 4
Symbol:	TMED4
UniProt ID:	<a href="#">Q7Z7H5</a>
Synonyms:	ERS25; GMP25iso; HNFLF; p24a3; p24alpha3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)