


## Product datasheet for **GA116673**

### Human SLC35F1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | SLC35F1 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member F1                                                                                                                                                                                                                         |
| Symbol:       | SLC35F1                                                                                                                                                                                                                                                                                                          |
| Components:   | GA116673G1, SLC35F1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA116673G2, SLC35F1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA116673G3, SLC35F1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q5T1Q4</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | C6orf169; dJ230I3.1                                                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)