


Product datasheet for **GA116678**

Human KCTD20 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCTD20 CRISPRa kit - CRISPR gene activation of human potassium channel tetramerization domain containing 20
Symbol:	KCTD20
UniProt ID:	Q7Z5Y7
Synonyms:	C6orf69; dj108K11.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)