

Product datasheet for **GA116682**

Human ZSCAN23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZSCAN23 CRISPRa kit - CRISPR gene activation of human zinc finger and SCAN domain containing 23
Symbol:	ZSCAN23
Components:	GA116682G1, ZSCAN23 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116682G2, ZSCAN23 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116682G3, ZSCAN23 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3MJ62
Synonyms:	dj29K1.3; dj29K1.3.1; ZNF390; ZNF453
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)