

Product datasheet for **GA116726**

Human PELI3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PELI3 CRISPRa kit - CRISPR gene activation of human pellino E3 ubiquitin protein ligase family member 3
Symbol:	PELI3
UniProt ID:	<u>Q8N2H9</u>
Synonyms:	MGC35521
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)