

Product datasheet for **GA116736**

Human HEPACAM2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HEPACAM2 CRISPRa kit - CRISPR gene activation of human HEPACAM family member 2
Symbol:	HEPACAM2
Components:	GA116736G1, HEPACAM2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116736G2, HEPACAM2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116736G3, HEPACAM2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A8MWW5
Synonyms:	MIKI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)