


Product datasheet for **GA116738**

Human C22ORF30 (PRR14L) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PRR14L CRISPRa kit - CRISPR gene activation of human proline rich 14 like
Symbol:	PRR14L
Components:	GA116738G1, C22ORF30 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116738G2, C22ORF30 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116738G3, C22ORF30 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5THK1</u>
Synonyms:	C22orf30
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)