

Product datasheet for **GA116744**

Human IPMK activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | IPMK CRISPRa kit - CRISPR gene activation of human inositol polyphosphate multikinase |
| Symbol: | IPMK |
| UniProt ID: | <u>Q8NFU5</u> |
| Synonyms: | Impk; inositol polyphosphate kinase 2; inositol polyphosphate multikinase; Ipk2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)