

Product datasheet for **GA116763**

Human POLDIP1 (KCTD13) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCTD13 CRISPRa kit - CRISPR gene activation of human potassium channel tetramerization domain containing 13
Symbol:	KCTD13
UniProt ID:	Q8WZ19
Synonyms:	BACURD1; FKSG86; hBACURD1; PDIP1; POLDIP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)