

Product datasheet for **GA116768**

Human **LRRC43** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | LRRC43 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 43 |
| Symbol: | LRRC43 |
| UniProt ID: | <u>Q8N309</u> |
| Synonyms: | MGC35140 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)