

Product datasheet for **GA116775**

Human TSGA10IP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSGA10IP CRISPRa kit - CRISPR gene activation of human testis specific 10 interacting protein
Symbol:	TSGA10IP
Components:	GA116775G1, TSGA10IP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116775G2, TSGA10IP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116775G3, TSGA10IP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3SY00
Synonyms:	FAM161C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)