

Product datasheet for **GA116805**

Human MPV17L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MPV17L CRISPRa kit - CRISPR gene activation of human MPV17 mitochondrial inner membrane protein like
Symbol:	MPV17L
Components:	GA116805G1, MPV17L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116805G2, MPV17L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116805G3, MPV17L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q2QL34
Synonyms:	M-LPH; MLPH1; MLPH2; MPV17L1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)