

Product datasheet for **GA116807**

Human DOS (CBARP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CBARP CRISPRa kit - CRISPR gene activation of human CACN subunit beta associated regulatory protein
Symbol:	CBARP
Components:	GA116807G1, DOS gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116807G2, DOS gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116807G3, DOS gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N350</u>
Synonyms:	BARP; C19orf26; DOS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)