

Product datasheet for **GA116808**

Human TAC4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAC4 CRISPRa kit - CRISPR gene activation of human tachykinin 4
Symbol:	TAC4
Components:	GA116808G1, TAC4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116808G2, TAC4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116808G3, TAC4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86UU9</u>
Synonyms:	EK; HK-1; HK1; PPT-C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)