

Product datasheet for **GA116810**

Human TMCO4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMCO4 CRISPRa kit - CRISPR gene activation of human transmembrane and coiled-coil domains 4
Symbol:	TMCO4
Components:	GA116810G1, TMCO4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116810G2, TMCO4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116810G3, TMCO4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5TGY1</u>
Synonyms:	RP5-1056L3.6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)