

Product datasheet for **GA116816**

Human **LRRC57** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	LRRC57 CRISPRa kit - CRISPR gene activation of human leucine rich repeat containing 57
Symbol:	LRRC57
Components:	GA116816G1, LRRC57 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116816G2, LRRC57 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116816G3, LRRC57 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N9N7</u>
Synonyms:	DKFZp686H1865; FLJ36812
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)