

Product datasheet for **GA116821**

Human Transmembrane protein 211 (TMEM211) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM211 CRISPRa kit - CRISPR gene activation of human transmembrane protein 211
Symbol:	TMEM211
UniProt ID:	Q6ICI0
Synonyms:	bA9F11.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)