

Product datasheet for **GA116830**

Human COL24A1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	COL24A1 CRISPRa kit - CRISPR gene activation of human collagen type XXIV alpha 1 chain
Symbol:	COL24A1
Components:	GA116830G1, COL24A1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116830G2, COL24A1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116830G3, COL24A1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q17RW2</u>
Synonyms:	MGC142214
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)