

Product datasheet for **GA116832**

Human PCSK9 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PCSK9 CRISPRa kit - CRISPR gene activation of human proprotein convertase subtilisin/kexin type 9
Symbol:	PCSK9
UniProt ID:	Q8NBP7
Synonyms:	FH3; HCHOLA3; LDLCQ1; NARC-1; NARC1; PC9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)