

## Product datasheet for **GA116837**

### Human C3orf43 (SMCO1) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SMCO1 CRISPRa kit - CRISPR gene activation of human single-pass membrane protein with coiled-coil domains 1
Symbol:	SMCO1
Components:	GA116837G1, C3orf43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116837G2, C3orf43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116837G3, C3orf43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q147U7</u>
Synonyms:	C3orf43
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)