

## Product datasheet for **GA116842**

### Human **PAN3** activation kit by **CRISPRa**

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                         |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                            |
| <b>Description:</b>  | PAN3 CRISPRa kit - CRISPR gene activation of human poly(A) specific ribonuclease subunit PAN3                                                                                                                                                                                                           |
| <b>Symbol:</b>       | PAN3                                                                                                                                                                                                                                                                                                    |
| <b>Components:</b>   | GA116842G1, PAN3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA116842G2, PAN3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA116842G3, PAN3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>UniProt ID:</b>   | <u><a href="#">Q58A45</a></u>                                                                                                                                                                                                                                                                           |
| <b>Synonyms:</b>     | OTTHUMP00000018179; OTTHUMP00000018180; OTTHUMP00000042400; PABP-dependent pol; PABP1-dependent poly A-specific ribonuclease subunit PAN3; PAN3 poly(A) specific ribonuclease subunit homolog (S. cerevisiae); PAN3 polyA specific ribonuclease subunit homolog                                         |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                 |



