

Product datasheet for **GA116859**

Human EML3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EML3 CRISPRa kit - CRISPR gene activation of human EMAP like 3
Symbol:	EML3
Components:	GA116859G1, EML3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116859G2, EML3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116859G3, EML3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q32P44</u>
Synonyms:	ELP95
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)