

Product datasheet for **GA116865**

Human TMEM151A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM151A CRISPRa kit - CRISPR gene activation of human transmembrane protein 151A
Symbol:	TMEM151A
Components:	GA116865G1, TMEM151A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116865G2, TMEM151A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116865G3, TMEM151A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8N4L1</u>
Synonyms:	TMEM151
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)