

Product datasheet for **GA116877**

Human KANK3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KANK3 CRISPRa kit - CRISPR gene activation of human KN motif and ankyrin repeat domains 3
Symbol:	KANK3
Components:	GA116877G1, KANK3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116877G2, KANK3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116877G3, KANK3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NY19</u>
Synonyms:	ANKRD47
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)