

Product datasheet for **GA116879**

Human SUNC1 (SUN3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SUN3 CRISPRa kit - CRISPR gene activation of human Sad1 and UNC84 domain containing 3
Symbol:	SUNC1
Components:	GA116879G1, SUNC1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116879G2, SUNC1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116879G3, SUNC1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8TAQ9</u>
Synonyms:	SUNC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)