

Product datasheet for **GA116893**

Human SHPRH activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SHPRH CRISPRa kit - CRISPR gene activation of human SNF2 histone linker PHD RING helicase
Symbol:	SHPRH
UniProt ID:	<u>Q149N8</u>
Synonyms:	bA54515.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)