


Product datasheet for **GA116897**

Human SNX33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SNX33 CRISPRa kit - CRISPR gene activation of human sorting nexin 33
Symbol:	SNX33
Components:	GA116897G1, SNX33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA116897G2, SNX33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA116897G3, SNX33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WV41
Synonyms:	SH3PX3; SH3PXD3C; SNX30
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)