

Product datasheet for **GA116917**

Human **TAS2R40** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TAS2R40 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 40
Symbol:	TAS2R40
UniProt ID:	P59535
Synonyms:	GPR60; T2R40; T2R58
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)