

Product datasheet for **GA116918**

Human **TAS2R41** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TAS2R41 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 41
Symbol:	TAS2R41
UniProt ID:	P59536
Synonyms:	T2R41; T2R59
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)