


## Product datasheet for **GA116919**

### Human **TAS2R43** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	TAS2R43 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 43
Symbol:	TAS2R43
UniProt ID:	<a href="#">P59537</a>
Synonyms:	T2R43; T2R52
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)