

Product datasheet for **GA116920**

Human **TAS2R44 (TAS2R31)** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAS2R31 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 31
Symbol:	TAS2R31
UniProt ID:	P59538
Synonyms:	T2R31; T2R44; T2R53; TAS2R44
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

