

Product datasheet for **GA116922**

Human **TAS2R46** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | TAS2R46 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 46 |
| Symbol: | TAS2R46 |
| UniProt ID: | P59540 |
| Synonyms: | T2R46; T2R54 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)