

## Product datasheet for **GA116925**

### Human **TAS2R49 (TAS2R20)** activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	TAS2R20 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 20
Symbol:	TAS2R20
UniProt ID:	<a href="#">P59543</a>
Synonyms:	T2R20; T2R49; T2R56; TAS2R49
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)