

Product datasheet for **GA116926**

Human T2R50 (TAS2R50) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TAS2R50 CRISPRa kit - CRISPR gene activation of human taste 2 receptor member 50
Symbol:	TAS2R50
UniProt ID:	P59544
Synonyms:	T2R50; T2R51; TAS2R51
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)