

Product datasheet for **GA116931**

Human **MAGI3** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | MAGI3 CRISPRa kit - CRISPR gene activation of human membrane associated guanylate kinase, WW and PDZ domain containing 3 |
| Symbol: | MAGI3 |
| UniProt ID: | Q5TCQ9 |
| Synonyms: | dj730K3.2; MAGI-3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)